

NISUS[®] Negative Pressure Wound Therapy System

A powerful yet portable device, allowing patient to achieve results while still maintaining mobility.

- ✓ FDA approved for use in hospitals, long-term care facilities, and home-care settings
- ✓ Lightweight and portable, with up to 24 hrs of battery life (can extend to 48 hrs with sleep mode)
- ✓ User-friendly interactive interface with easy-to-read LCD screen



NPWT Wound Dressing Kit

Nisus NPWT Pump

Nisus NPWT Canister(s)

The Cork Medical Products Nisus Negative Pressure Wound Therapy System is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

NISUS NPWT SYSTEM FEATURES

The Cork Medical Products Nisus Negative Pressure Wound Therapy System is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- **Recommended for patients with:**
 - Chronic, acute, traumatic, sub-acute and dehiscent wounds
 - Partial-thickness burns
 - Pressure ulcers, diabetic ulcers, and venous ulcers
 - Flaps and grafts
- **Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage**
- **IV pole accessible**
- **Powerful flow rate of ~10 L/min**
- **Multiple unique features for clinical enhancement of negative pressure wound therapy such as:**
 - Patient runtime
 - Sleep feature
 - Alarm log
 - Lock restriction
 - Variable Intermittent and Continuous Therapy Modes
- **Convenient rear compartment battery access and replacement**

COMPONENT SPECS

Comprised of three essential components that actively work together to promote wound healing, the Nisus NPWT System includes the Nisus NPWT Pump, the Nisus NPWT Canisters, and the Cork NPWT Wound Dressing Kit.

Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical Products components with the Nisus Negative Pressure Wound Therapy System.

- **Battery charger and neoprene carrying bag included**
- **Medicare HCPCS Code: E2402**

Model No.	Component
CMPP-100	Nisus NPWT Pump
CPC-250	Nisus NPWT 250mL Canister
CPC-500	Nisus NPWT 500mL Canister
WDK-250	NPWT Wound Dressing Kit (Small)
WDK-300	NPWT Wound Dressing Kit (Medium)
WDK-400	NPWT Wound Dressing Kit (Large)
CLEANSE-KIT-A-18	NPWT Cleanse Therapy Kit
NIS-BRDR-A-16	Bridge Foam
FOAM-WHITE-A-16	HydroLite White Foam
CSP-100-ES-EXT	24" Extension Tubing
CWY-100	Y-Connector



In Partnership with

CORK
MEDICAL

**For more information, contact the InfuSystem NPWT Team at
1-844-534-1997 or npwt@infusystem.com**

InfuSystem[®]
SAFE. SMART. TRUSTED.SM