

Nisus[®] Negative Pressure Wound Therapy (NPWT) Canisters

The Nisus NPWT Canisters are designed to be used with Cork Medical Products Nisus Negative Pressure Wound Therapy Pump to provide negative pressure wound therapy to aid and promote wound healing by removing excess exudates, infectious material and tissue debris.



- ✓ Latch design provides streamlined, one-handed canister release and removal
- ✓ Graduated indications on canister aid in measuring exudate fluid

NISUS NPWT CANISTERS PRODUCT FEATURES

The Nisus NPWT Canisters, used in conjunction with the Nisus Negative Pressure Wound Therapy Pump and the Wound Dressing Kit, are indicated for use in patients who would benefit from negative pressure wound therapy. The Nisus unique canister design sustains pressure under variable conditions, allowing patients the flexibility to be portable while maintaining a hermetic seal.

- **Canister latch design** features a hermetic seal that can sustain pressure under variable conditions
- **User-friendly design** provides streamlined, one-handed canister release and removal
- **State-of-the-art dam design** keeps membrane dry and prevents false canister full alarm
- **Graduated indications** on canister aid in measuring exudate fluid
- **Available in two (2) volume sizes:** 250 mL for portability and 500 mL for larger wound drainage
- **Canister is preassembled** with a 30" drainage tube with clamp and male luer for easy connection

COMPONENT SPECIFICATIONS

The Nisus NPWT Canisters work exclusively with the Cork Medical Products NPWT Wound Therapy Pump and Dressing Kit to provide a safe complete wound care system.

- Single-use and housed in a non-sterile Tyvek peel pouch, canisters collect and store wound drainage
- Medicare HCPCS Code: A7000

 **500 mL Canister for larger wounds**
250 mL Canister for portability



Component	Model No.
Nisus NPWT Pump	CMPP-100
Nisus NPWT Canisters	CPC-250 (250 mL) or CPC-500 (500 mL)
NPWT Wound Dressing Kit	WDK-250 (Small), WDK-300 (Medium), or WDK-400 (Large)

Please use under the advisement of a physician and carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical Products components with the Nisus Negative Pressure Wound Therapy Canisters.

 In Partnership with


**For more information, contact the InfuSystem NPWT Team at
1-844-534-1997 or npwt@infusystem.com**

